

β-Catenin Antibody

Catalog No: #47993

Package Size: #47993-1 50ul #47993-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	β-Catenin Antibody
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB, ICC, IHC
Species Reactivity	Hu,Ms,Rt, zebrafish
Immunogen Description	peptide
Conjugates	Unconjugated
Other Names	Beta catenin antibody Beta-catenin antibody Cadherin associated protein antibody Catenin (cadherin associated protein), beta 1, 88kDa antibody Catenin beta 1 antibody Catenin beta-1 antibody CATNB antibody CHBCAT antibody CTNB1_HUMAN antibody CTNNB antibody CTNNB1 antibody DKFZp686D02253 antibody FLJ25606 antibody FLJ37923 antibody OTTHUMP00000162082 antibody OTTHUMP00000165222 antibody OTTHUMP00000165223 antibody OTTHUMP00000209288 antibody OTTHUMP00000209289 antibody
Accession No.	Swiss-Prot#:Q02248
Calculated MW	86 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details


WB 1:1,000

ICC: 1:200

Images



Western blot analysis on cell lysates using anti-β-Catenin polyclonal antibody



ICC staining β -Catenin (green) in HeLa cells. Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Background

β -Catenin is a multifunctional adaptor protein involved in cadherin-mediated cell-cell adhesion and in responding to the activation of several signal transduction pathways. CK1 phosphorylates β -catenin on Ser45. This phosphorylation event primes β -catenin for subsequent phosphorylation by GSK-3. GSK-3 β destabilizes β -catenin by phosphorylating it at Ser33, Ser37 and Thr41.). Mutations in β -Catenin have been implicated in its accumulation and various forms of carcinomas.

References

Published Papers

el at., Circular RNA hsa_circ_0003823 promotes the Tumor Progression, Metastasis and Apatinib Resistance of Esophageal Squamous Cell Carcinoma by miR-607/CRISP3 Axis. In Int J Biol Sci on 2022 Sep 21 by Yu-Ming Wang, Qi-Wu Zhao, et al..PMID:36263172, , (2022)
[PMID:36263172](#)

Zhu Tian-Yi;Hu Po;Mi Yu-Hui;Zhang Jun-Li;Xu An-Na;Gao Ming-Tong;Zhang Ying-Ying;Shen San-Bing;Yang Guang-Ming;Pan Yang el at., Telomerase reverse transcriptase gene knock-in unleashes enhanced longevity and accelerated damage repair in mice, , (2024)
[PMID:](#)

Zhipeng Zhang;Yanmei Wu;Wenqiang Liang;Zhifang Liao;Hongbo Liao;Xingxing Xing;Wenxin Yi;Zixuan Liu;Yicheng Li;Mengya Shi;Dongling Lin;Ting Gu;Biao Wu;Mingzhi Zou;Huilai Miao;Xin Wu el at., Eurycomalactone switched hepatocellular carcinoma cells into quiescence through 5'tRFAla/DVL/ β -catenin pathway inhibition., , (2025)
[PMID:40128187](#)

Jie Guo;Suhan Zhou;Honghong Wang;Xingyu Qiu;Fang Dong;Shan Jiang;Nan Xu;Yu Cui;Ruisheng Liu;Pengyun Li;Zufu Ma;Liang Zhao;En Yin Lai el at., ADAMTS13 attenuates renal fibrosis by suppressing thrombospondin 1 mediated TGF- β 1/Smad3 activation., , (2025)
[PMID:39929281](#)

Shuting Gao;Huihua Li;Zekun Li;Hong Wang;Xinyue Li;Shengyan Yang;Lin Huang;Baoping Zhang;Kailiang Zhang;James Kit Hon Tsoi;Jian He;Waruna Lakmal Dissanayaka el at., Multifunctional Injectable Bioadhesive with Toll-like Receptor 4 and Myeloid Differentiation Factor 2 Antagonistic Anti-inflammatory Potential for Periodontal Regeneration., , (2025)
[PMID:39951685](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.