

Human IgG heavy chain Antibody

Catalog No: #48056

Package Size: #48056-1 50ul #48056-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Human IgG heavy chain Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purified
Applications	WB
Species Reactivity	Hu
Immunogen Description	Human IgG
Conjugates	Unconjugated
Accession No.	Swiss-Prot#:
Calculated MW	50 kDa
Formulation	1*TBS (pH7.4), 0.5%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:5,000

Images



Western blot analysis on human IgG protein using anti-Human IgG heavy chain Mouse mAb (Cat. # M0806-1).

References

Note: This product is for in vitro research use only and is not intended for use in humans or animals.