

Rabbit IgG κ Light Chain Antibody

Catalog No: #48059

Package Size: #48059-1 50ul #48059-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Rabbit IgG κ Light Chain Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purified
Applications	WB, IP
Species Reactivity	Rabbit
Immunogen Description	Rabbit IgG
Conjugates	Unconjugated
Accession No.	Swiss-Prot#:
Calculated MW	25 kDa
Formulation	1*TBS (pH7.4), 0.5%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

IP: 1:500-1,000

Images



Western blot analysis on rabbit IgG protein using anti-Rabbit IgG kappa light chain Mouse mAb (Cat. # M0809-1).

References

Note: This product is for in vitro research use only and is not intended for use in humans or animals.