

Histone Deacetylase 3(HDAC3) Antibody

Catalog No: #48123

Package Size: #48123-1 50ul #48123-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone Deacetylase 3(HDAC3) Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	B2-G9
Purification	Peptide affinity purified
Applications	WB, ICC, IHC
Species Reactivity	Hu, Ms, Rt
Immunogen Description	peptide
Conjugates	Unconjugated
Accession No.	Swiss-Prot#:O15379
Calculated MW	48 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

References

Note: This product is for in vitro research use only and is not intended for use in humans or animals.