

VSV-g Tag Rabbit mAb

Catalog No: #48191

Package Size: #48191-1 50ul #48191-2 100ul


Description

Product Name	VSV-g Tag Rabbit mAb
Clonality	Monoclonal
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Recombinant protein
Immunogen Type	Protein
Immunogen Description	Recombinant Protein of VSV G tag
Conjugates	Unconjugated
Other Names	Vesicular stomatitis virus glycoprotein tag antibody VSV epitope tag antibody VSV tag antibody VSV-G epitope tag antibody VSVG antibody vsvg tag antibody
Accession No.	Swiss-Prot#:
Calculated MW	Recombinant protein dependent
Concentration	0.3 mg/ml
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Application Details

WB 1:1000-1:10000

Images



Western blot analysis of extracts from VSV-G tag fusion protein, using VSV-G tag antibody.

Background

Vesicular stomatitis virus (VSV), an enveloped RNA virus from the Rhabdoviridae family, is released from the plasma membrane of host cells by a process called budding. The glycoprotein (VSV-g) contains a domain in its extracellular membrane proximal stem that appears to be needed for efficient VSV budding.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.