

PLC β 2 Polyclonal Antibody

Catalog No: #48323

Package Size: #48323-1 50ul #48323-2 100ul


Description

Product Name	PLC β 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC
Species Reactivity	Human;Mouse;Rat
Specificity	PLC β 2 Polyclonal Antibody detects endogenous levels of PLC β 2 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human PLC β 2.
Conjugates	Unconjugated
Target Name	PLCB2
Other Names	PLCB2; 1-phosphatidylinositol 4; 5-bisphosphate phosphodiesterase beta-2; Phosphoinositide phospholipase C-beta-2; Phospholipase C-beta-2; PLC-beta-2
Accession No.	Q00722
Calculated MW	150kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000;IHC-p 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung.

Published Papers

el at., Salicin alleviates periodontitis via Tas2r143/gustducin signaling in fibroblasts. In Front Immunol on 2024 Mar 28 by Zhiying Zhang, Zhiyan Zhou,et al..PMID:38605968, , (2024)

[PMID:38605968](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.