

PD1 Blocking Antibody

Catalog No: #48361

Package Size: #48361-1 50ul #48361-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PD1 Blocking Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3A11
Purification	ProG affinity purified
Applications	FC
Species Reactivity	Hu, Ms
Immunogen Description	Recombinant protein
Conjugates	Unconjugated
Accession No.	Swiss-Prot#:Q15116
Calculated MW	32 kDa
Formulation	1*TBS (pH7.4), 0.5%BSA, 50%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

FC: 1:50-1:100

Images



Flow cytometric analysis of Jurkat cells with PD1 blocking antibody at 1/100 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black).

References

Note: This product is for in vitro research use only and is not intended for use in humans or animals.