

KLHL11 Antibody

Catalog No: #48456

Package Size: #48456-1 50ul #48456-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KLHL11 Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	A8-B9
Purification	ProA affinity purified
Applications	WB, ICC/IF, FC
Species Reactivity	Hu
Immunogen Description	Recombinant protein
Conjugates	Unconjugated
Other Names	FLJ10572 antibody Kelch like 11 (Drosophila) antibody Kelch like 11 antibody Kelch like protein 11 antibody Kelch-like protein 11 antibody KLH11_HUMAN antibody KLHL 11 antibody Klh11 antibody
Accession No.	Swiss-Prot#:Q9NVR0
Calculated MW	80 kDa
Formulation	1*TBS (pH7.4), 1%BSA, Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:500-1:1,000 ICC: 1:50-1:200 FC: 1:50-1:200

Background

Component of a cullin-RING-based BCR (BTB-CUL3-RBX1) E3 ubiquitin-protein ligase complex that mediates the ubiquitination of target proteins, leading most often to their proteasomal degradation.

References

Note: This product is for in vitro research use only and is not intended for use in humans or animals.