

SLC25A20 rabbit pAb

Catalog No: #48787

Package Size: #48787-1 50ul #48787-2 100ul

Description

Product Name	SLC25A20 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of SLC25A20 at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human SLC25A20
Conjugates	Unconjugated
Target Name	SLC25A20
Calculated MW	33kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.507% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.