

TRPV3 Polyclonal Antibody

Catalog No: #48836

Package Size: #48836-1 50ul #48836-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TRPV3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse
Specificity	TRPV3 Polyclonal Antibody detects endogenous levels of TRPV3 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human TRPV3.
Conjugates	Unconjugated
Target Name	TRPV3
Other Names	TRPV3; Transient receptor potential cation channel subfamily V member 3; TrpV3; Vanilloid receptor-like 3; VRL-3
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC-p: 1:100-300 ELISA: 1:20000.. IF 1:50-200

Images



Western Blot analysis of extracts from CACO-2 cells, using TRPV3 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:100

Note: This product is for in vitro research use only and is not intended for use in humans or animals.