

T2R5 Polyclonal Antibody

Catalog No: #48849

Package Size: #48849-1 50ul #48849-2 100ul


Description

Product Name	T2R5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	T2R5 Polyclonal Antibody detects endogenous levels of T2R5 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human TAS2R5.
Conjugates	Unconjugated
Target Name	TAS2R5
Other Names	TAS2R5; Taste receptor type 2 member 5; T2R5
Calculated MW	35kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000. ELISA: 1:5000.

Images



Western blot analysis of the lysates from HepG2 cells using TAS2R5 antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.