

TIS11B Polyclonal Antibody

Catalog No: #48856

Package Size: #48856-1 50ul #48856-2 100ul


Description

Product Name	TIS11B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	TIS11B Polyclonal Antibody detects endogenous levels of TIS11B protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human TISB.
Conjugates	Unconjugated
Target Name	ZFP36L1
Other Names	ZFP36L1; BERG36; BRF1; ERF1; RNF162B; TIS11B; Zinc finger protein 36; C3H1 type-like 1; Butyrate response factor 1; EGF-response factor 1; ERF-1; Protein TIS11B
Calculated MW	36kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:10000.. IF 1:50-200

Images



Western Blot analysis of various cells using TIS11B Polyclonal Antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using TISB Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.