

MT1A rabbit pAb

Catalog No: #48877

Package Size: #48877-1 50ul #48877-2 100ul


Description

Product Name	MT1A rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF
Species Reactivity	Human;Rat;Mouse
Specificity	This antibody detects endogenous levels of MT1A at Human
Immunogen Description	Synthesized peptide derived from human MT1A
Conjugates	Unconjugated
Target Name	MT1A MT1S
Calculated MW	7kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC-p 1oO 50-200. IF 1:50-200

Images



Immunohistochemical analysis of paraffin-embedded human liver cancer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.