

## LCE1F rabbit pAb

Catalog No: #48879

Package Size: #48879-1 50ul #48879-2 100ul


## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | LCE1F rabbit pAb                                                                                                      |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Isotype               | IgG                                                                                                                   |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB                                                                                                                    |
| Species Reactivity    | Human;Rat;Mouse                                                                                                       |
| Specificity           | This antibody detects endogenous levels of LCE1F at Human                                                             |
| Immunogen Description | Synthesized peptide derived from human LCE1F                                                                          |
| Conjugates            | Unconjugated                                                                                                          |
| Target Name           | LCE1F LEP6                                                                                                            |
| Calculated MW         | 13kD                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C/1 year                                                                                                 |

## Application Details

WB 1αO 500-2000

## Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.