

TIS11D Polyclonal Antibody

Catalog No: #48901

Package Size: #48901-1 50ul #48901-2 100ul


Description

Product Name	TIS11D Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse
Specificity	TIS11D Polyclonal Antibody detects endogenous levels of TIS11D protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human TISD.
Conjugates	Unconjugated
Target Name	ZFP36L2
Other Names	ZFP36L2; BRF2; ERF2; RNF162C; TIS11D; Zinc finger protein 36; C3H1 type-like 2; ZFP36-like 2; Butyrate response factor 2; EGF-response factor 2; ERF-2; Protein TIS11D
Calculated MW	55kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:20000.. IF 1:50-200

Images



Western Blot analysis of various cells using TIS11D Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human small intestinal carcinoma tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.