

S2P Polyclonal Antibody

Catalog No: #48919

Package Size: #48919-1 50ul #48919-2 100ul


Description

Product Name	S2P Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB;IP;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	S2P Polyclonal Antibody detects endogenous levels of S2P protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human MBTPS2.
Conjugates	Unconjugated
Target Name	MBTPS2
Other Names	MBTPS2; S2P; Membrane-bound transcription factor site-2 protease; Endopeptidase S2P; Sterol regulatory element-binding proteins intramembrane protease; SREBPs intramembrane protease
Calculated MW	57kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000 IF 1:200 - 1:1000. ELISA: 1:20000

Images



Immunofluorescence analysis of A549 cells, using MBTPS2 Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.