

TRPC6 rabbit pAb

Catalog No: #48948

Package Size: #48948-1 50ul #48948-2 100ul


Description

Product Name	TRPC6 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of TRPC6 at Human/Mouse
Immunogen Description	Synthesized peptide derived from human TRPC6
Conjugates	Unconjugated
Target Name	TRPC6 TRP6
Calculated MW	102kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1αO 500-2000

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.