

GPHB5 Rabbit pAb

Catalog No: #48970

Package Size: #48970-1 50ul #48970-2 100ul



Description

Product Name	GPHB5 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	IHC-P, IHC-F, ICC, IF, ELISA
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human GPHB5: 25-100/130
Conjugates	Unconjugated
Other Names	CGH; Glycoprotein hormone beta 5; GPHB5_HUMAN; Glycoprotein hormone beta 5 precursor; Glycoprotein hormone beta subunit; GPB5; Thyrostimulin subunit beta; ZLUT1.
Calculated MW	12kD
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P 1:100-500
IHC-F 1:100-500
ICC 1:100-500
IF 1:100-500
ELISA 1:5000-10000

Published Papers

Shiding Li;Hao Sun;Fei Fang;Siyi Zhang;Junzhao Chen;Chunyi Shao;Yao Fu;Liangbo Chen et al., The Spatial Transcriptomic Atlas of Human Limbus and Vital Niche Microenvironment Regulating the Fate of Limbal Epithelial Stem Cells., , (2025)
[PMID:40131296](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.