

ALDH3A1 Rabbit pAb

Catalog No: #49026

Package Size: #49026-1 50ul #49026-2 100ul

Description

Product Name	ALDH3A1 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	WB, IHC-P, IHC-F, ICC, IF, ELISA
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human ALDH3A1: 151-250/453
Conjugates	Unconjugated
Other Names	AHD 4; AHD C; Ahd4; AHDC; AL3A1_HUMAN; Aldehyde dehydrogenase 3; Aldehyde dehydrogenase 3 family member A1; Aldehyde dehydrogenase 3 family, member A1; Aldehyde dehydrogenase 3A1; Aldehyde dehydrogenase; Aldehyde dehydrogenase class 3; Aldehyde dehydrogenase dimeric NADP preferring; Aldehyde dehydrogenase family 3 member A1; Aldehyde dehydrogenase family 3 subfamily A1; Aldehyde dehydrogenase isozyme 3; Aldehyde dehydrogenase type III; Aldehyde dehydrogenase, dimeric NADP preferring; ALDH 3; ALDH 3A1; ALDH III; ALDH3 A1; ALDH3; ALDH3A 1; Aldh3a1; ALDHIII; dimeric NADP-preferring; MGC10406; Stomach aldehyde dehydrogenase; Tumor associated aldehyde dehydrogenase tumor ALDH.
Calculated MW	50kD
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

IHC-P 1:100-500

IHC-F 1:100-500

ICC 1:100-500

IF 1:100-500

ELISA 1:5000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.