

ASCL1 Polyclonal Antibody

Catalog No: #49035

Package Size: #49035-1 50ul #49035-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ASCL1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	ASCL1 Polyclonal Antibody detects endogenous levels of ASCL1 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human ASCL1.
Conjugates	Unconjugated
Target Name	ASCL1
Other Names	ASCL1;ASH1;BHLHA46;HASH1;Achaete-scute homolog 1;ASH-1;hASH1;Class A basic helix-loop-helix protein 46;bHLHa46
Calculated MW	25kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:40000.. IF 1:50-200

Images



Western blot analysis of lysate from Jurkat cells, using ASCL1 antibody.



Immunohistochemistry analysis of ASCL1 antibody in paraffin-embedded human brain tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.