

MIF Polyclonal Antibody

Catalog No: #49061

Package Size: #49061-1 50ul #49061-2 100ul


Description

Product Name	MIF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IF;WB;IHC;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	MIF Polyclonal Antibody detects endogenous levels of MIF protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human MIF.
Conjugates	Unconjugated
Target Name	MIF
Other Names	MIF;GLIF;MMIF;Macrophage migration inhibitory factor;MIF;Glycosylation-inhibiting factor;GIF;L-dopachrome isomerase;L-dopachrome tautomerase;Phenylpyruvate tautomerase
Calculated MW	12kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

IF 1:50-200 WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300


Images



Western blot analysis of lysate from HepG2 cells, using MIF antibody.



Immunofluorescence analysis of A549.



Immunohistochemical analysis of paraffin-embedded Mouse-kidney tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.