

## DNAJB9 Rabbit pAb

Catalog No: #49066

Package Size: #49066-1 50ul #49066-2 100ul

## Description

|                       |                                                                                                                                                                                                                                                                                                                               |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | DNAJB9 Rabbit pAb                                                                                                                                                                                                                                                                                                             |
| Host Species          | Rabbit                                                                                                                                                                                                                                                                                                                        |
| Clonality             | Polyclonal                                                                                                                                                                                                                                                                                                                    |
| Isotype               | IgG                                                                                                                                                                                                                                                                                                                           |
| Purification          | affinity purified by Protein A                                                                                                                                                                                                                                                                                                |
| Applications          | WB, IHC-P, IHC-F, ICC, IF, ELISA                                                                                                                                                                                                                                                                                              |
| Species Reactivity    | Human;Mouse;Rat                                                                                                                                                                                                                                                                                                               |
| Immunogen Description | KLH conjugated synthetic peptide derived from human DNAJB9: 24-120/223                                                                                                                                                                                                                                                        |
| Conjugates            | Unconjugated                                                                                                                                                                                                                                                                                                                  |
| Other Names           | DKFZP564F1862; DnaJ (Hsp40) homolog subfamily B member 9; DnaJ homolog subfamily B member 9; DNAJB9; DNJB9_HUMAN; Endoplasmic reticulum DnaJ homolog 4; ERdj4; MDG 1; Mdg-1; Microvascular endothelial differentiation gene 1; Microvascular endothelial differentiation gene 1 protein; MST049; MSTP049; OTTHUMP00000206840. |
| Calculated MW         | 23 kDa                                                                                                                                                                                                                                                                                                                        |
| Concentration         | 1mg/ml                                                                                                                                                                                                                                                                                                                        |
| Formulation           | 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.                                                                                                                                                                                                                                                             |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                                                                                                                                                                                                                               |

## Application Details

WB 1:500-2000

IHC-P 1:100-500

IHC-F 1:100-500

ICC 1:100-500

IF 1:100-500

ELISA 1:5000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.