

SCNM1 Polyclonal Antibody

Catalog No: #49083

Package Size: #49083-1 50ul #49083-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SCNM1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	SCNM1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Conjugates	Unconjugated
Target Name	SCNM1
Calculated MW	25kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000 ELISA 1:5000-20000

Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.