

MUC6 rabbit pAb

Catalog No: #49094

Package Size: #49094-1 50ul #49094-2 100ul


Description

Product Name	MUC6 rabbit pAb
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of Human MUC6
Immunogen Description	Synthesized peptide derived from human MUC6
Conjugates	Unconjugated
Target Name	MUC6
Other Names	Mucin-6 (MUC-6;Gastric mucin-6)
Calculated MW	268kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000;IHC-p 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human liver cancer

Note: This product is for in vitro research use only and is not intended for use in humans or animals.