

CR3L1 Polyclonal Antibody

Catalog No: #49153

Package Size: #49153-1 50ul #49153-2 100ul

Description

Product Name	CR3L1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	CR3L1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Conjugates	Unconjugated
Target Name	CREB3L1 OASIS PSEC0238
Calculated MW	57kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000 ELISA 1:5000-20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.