

METRL rabbit pAb

Catalog No: #49186

Package Size: #49186-1 50ul #49186-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	METRL rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of METRL at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human METRL
Conjugates	Unconjugated
Target Name	METRNL
Calculated MW	34kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-2000;IHC 1:50-300

Images



Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4°over night



Immunohistochemical analysis of paraffin-embedded human liver cancer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.