

## USP32 Polyclonal Antibody

Catalog No: #49189

Package Size: #49189-1 50ul #49189-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                                                                                                       |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | USP32 Polyclonal Antibody                                                                                                                                                             |
| Host Species          | Rabbit                                                                                                                                                                                |
| Clonality             | Polyclonal                                                                                                                                                                            |
| Isotype               | IgG                                                                                                                                                                                   |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                 |
| Applications          | IHC;IF;WB;ELISA                                                                                                                                                                       |
| Species Reactivity    | Human;Mouse                                                                                                                                                                           |
| Specificity           | USP32 Polyclonal Antibody detects endogenous levels of USP32 protein.                                                                                                                 |
| Immunogen Description | The antiserum was produced against synthesized peptide derived from human USP32.                                                                                                      |
| Conjugates            | Unconjugated                                                                                                                                                                          |
| Target Name           | USP32                                                                                                                                                                                 |
| Other Names           | USP32;USP10;Ubiquitin carboxyl-terminal hydrolase 32;Deubiquitinating enzyme 32;Renal carcinoma antigen NY-REN-60;Ubiquitin thioesterase 32;Ubiquitin-specific-processing protease 32 |
| Concentration         | 1 mg/ml                                                                                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                               |
| Storage               | Store at -20°C/1 year                                                                                                                                                                 |

## Application Details

WB 1:500-2000 IHC 1:100 - 1:300. ELISA: 1:40000.. IF 1:50-200

## Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using USP32 Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.