

TBX2 Polyclonal Antibody

Catalog No: #49236

Package Size: #49236-1 50ul #49236-2 100ul


Description

Product Name	TBX2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse
Specificity	TBX2 Polyclonal Antibody detects endogenous levels of TBX2 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human TBX2.
Conjugates	Unconjugated
Target Name	TBX2
Other Names	TBX2;T-box transcription factor TBX2;T-box protein 2
Calculated MW	74kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:20000.. IF 1:50-200

Images



Western blot analysis of lysates from HeLa cells, using TBX2 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human brain tumor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.