

## VPS52 rabbit pAb

Catalog No: #49241

Package Size: #49241-1 50ul #49241-2 100ul

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | VPS52 rabbit pAb                                                                                                      |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Isotype               | IgG                                                                                                                   |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB                                                                                                                    |
| Species Reactivity    | Human;Mouse;Rat                                                                                                       |
| Specificity           | This antibody detects endogenous levels of VPS52 at Human, Mouse,Rat                                                  |
| Immunogen Description | Synthesized peptide derived from human VPS52                                                                          |
| Conjugates            | Unconjugated                                                                                                          |
| Target Name           | VPS52 SACM2L                                                                                                          |
| Calculated MW         | 80kD                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C/1 year                                                                                                 |

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.