

CHAC1 Polyclonal Antibody

Catalog No: #49244

Package Size: #49244-1 50ul #49244-2 100ul


Description

| | |
|-----------------------|---|
| Product Name | CHAC1 Polyclonal Antibody |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | CHAC1 Polyclonal Antibody detects endogenous levels of protein. |
| Applications | WB;ELISA |
| Species Reactivity | Human;Rat;Mouse |
| Specificity | CHAC1 Polyclonal Antibody detects endogenous levels of protein. |
| Immunogen Description | Synthesized peptide derived from part region of human protein |
| Conjugates | Unconjugated |
| Target Name | CHAC1 BOTCH |
| Calculated MW | 24kD |
| Concentration | 1 mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

WB 1:500-2000 ELISA 1:5000-20000

Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.