

p14ARF Polyclonal Antibody

Catalog No: #49256

Package Size: #49256-1 50ul #49256-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	p14ARF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	p14 Polyclonal Antibody detects endogenous levels of p14 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human p14 ARF.
Conjugates	Unconjugated
Target Name	CDKN2A
Other Names	CDKN2A;CDKN2;MLM;Cyclin-dependent kinase inhibitor 2A; isoform 4;p14ARF;p19ARF
Calculated MW	18kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. IF 1:200 - 1:1000. ELISA: 1:20000.

Images



Western Blot analysis of various cells using p14 Polyclonal Antibody diluted at 1:500



Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using p14 ARF Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.