

FOXO4 Rabbit mAb

Catalog No: #49279

Package Size: #49279-1 50ul #49279-2 100ul


Description

Product Name	FOXO4 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Clone No.	JJ09-11
Purification	ProA affinity purified
Applications	WB, ICC/IF, IP, FC
Species Reactivity	Hu, Ms, Rt
Immunogen Description	recombinant protein
Conjugates	Unconjugated
Other Names	AFX antibody AFX1 antibody Afxh antibody ALL1-fused gene from X chromosome antibody Fork head domain transcription factor AFX1 antibody Forkhead box O4 antibody Forkhead box protein O4 antibody FOXO 4 antibody Foxo4 antibody FOXO4_HUMAN antibody MGC117660 antibody MGC120490 antibody Mixed lineage leukemia, translocated to, 7 antibody MLLT7 antibody Myeloid/lymphoid or mixed lineage leukemia (trithorax homolog, Drosophila); translocated to, 7 antibody Myeloid/lymphoid or mixed lineage leukemia, translocated to, 7 antibody RGD1561201 antibody
Accession No.	Swiss-Prot#:P98177
Calculated MW	54 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C


Application Details

WB: 1:1,000 ICC: 1:100-1:500 FC: 1:50-1:100


Images




ICC staining FOXO4 in 293T cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining FOXO4 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining FOXO4 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Flow cytometric analysis of 293 cells with FOXO4 antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody

Background

FKHR (for forkhead in rhabdomyosarcoma), FKHL1, and AFX1 are members of a subfamily of the forkhead family of transcription factors. AFX1, also known as FoxO4, is expressed in a wide variety of tissues and, like other FKHR proteins, AFX1 contains a single forkhead domain and serine-proline-rich region, which mediate DNA binding. AFX1-mediated transcriptional activation is regulated by the serine/threonine kinase Akt1, which phosphorylates AFX1 and in turn, sequesters AFX1 in the cytosol, thereby blocking nuclear localization and DNA binding. Genetic mutations in FKHR genes, including the t(2;13) and t(1;3) translocations, are commonly found in alveolar rhabdomyosarcomas. Additionally, the t(x;11) translocation of the AFX1 gene, which involves the fusion of a serine-proline-rich sequence of AFX1 to the carboxy terminus of a truncated MLL, results in acute lymphocytic leukemia.

References

Note: This product is for in vitro research use only and is not intended for use in humans or animals.