

LY6A Rabbit mAb

Catalog No: #49480

Package Size: #49480-1 50ul #49480-2 100ul



**SAB**<sup>®</sup>

Signalway Antibody

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

|                       |                                                                                                                                                                                                                                                                                                                    |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | LY6A Rabbit mAb                                                                                                                                                                                                                                                                                                    |
| Host Species          | Recombinant Rabbit                                                                                                                                                                                                                                                                                                 |
| Clonality             | Monoclonal                                                                                                                                                                                                                                                                                                         |
| Clone No.             | JM53-11                                                                                                                                                                                                                                                                                                            |
| Purification          | ProA affinity purified                                                                                                                                                                                                                                                                                             |
| Applications          | WB, IHC                                                                                                                                                                                                                                                                                                            |
| Species Reactivity    | Ms                                                                                                                                                                                                                                                                                                                 |
| Immunogen Description | recombinant protein                                                                                                                                                                                                                                                                                                |
| Conjugates            | Unconjugated                                                                                                                                                                                                                                                                                                       |
| Other Names           | Ly-6A.2/Ly-6E. antibody Ly-6A.2/Ly-6E.1 antibody Ly6a antibody LY6A_MOUSE antibody Ly6al antibody Lymphocyte antigen 6A-2/6E-1 antibody Lymphocyte antigen Ly 6A.2/Ly 6E.1 antibody SCA-1 antibody Stem cell antigen 1 antibody T cell activating protein antibody T-cell-activating protein antibody TAP antibody |
| Accession No.         | Swiss-Prot#:P05533                                                                                                                                                                                                                                                                                                 |
| Calculated MW         | 14/40 kDa                                                                                                                                                                                                                                                                                                          |
| Formulation           | 1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.                                                                                                                                                                                                                                               |
| Storage               | Store at -20°C                                                                                                                                                                                                                                                                                                     |


Application Details

WB: 1:1,000-5,000IHC: 1:50-1:200


Images



Western blot analysis of LY6A on mouse brain cells lysates using anti-LY6A antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded brain mouse tissue using anti-LY6A antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue using anti-LY6A antibody. Counter stained with hematoxylin.

## References

Note: This product is for in vitro research use only and is not intended for use in humans or animals.