

CLK2 Polyclonal Antibody

Catalog No: #49996

Package Size: #49996-1 50ul #49996-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLK2 Polyclonal Antibody
Isotype	IgG
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	CLK2 Polyclonal Antibody detects endogenous levels of CLK2 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human CLK2.
Conjugates	Unconjugated
Target Name	CLK2
Other Names	CLK2;Dual specificity protein kinase CLK2;CDC-like kinase 2
Calculated MW	60kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:20000.. IF 1:50-200

Images



Western blot analysis of lysates from COLO205 and A549 cells, using CLK2 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded Human brain.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.