

HDAC5(Phospho-Ser498) antibody blocking peptide

Catalog No: #51193

Package Size: #51193 100ul

Description

Product Name	HDAC5(Phospho-Ser498) antibody blocking peptide
Applications	NC
Conjugates	Unconjugated
Target Name	HDAC5
Modification	Phospho
Concentration	1.0mg/ml

Published Papers

Sheng Xia, Xiaogang Li, Teri Johnson et al., Polycystin-dependent fluid flow sensing targets histone deacetylase 5 to prevent the development of renal cysts., Development, 137, 1075-1084(2010)

[PMID:20181743](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.