

ABCA1 Mouse mAb

Catalog No: #52906

Package Size: #52906-1 50ul #52906-2 100ul



Description

Product Name	ABCA1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Purification	Affinity Purified
Applications	ICC/IF,FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ABCA1 expressed in E. Coli.
Conjugates	Unconjugated
Target Name	ABCA1
Other Names	TGD; ABC1; CERP; ABC-1; HDLDT1; HPALP1; HDLCQTL13
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IF: 1/50-1/200 FC: 1/50-1/100

Published Papers

Hao Xu;Hao Xu;Xueni Sun;Xueni Sun;Miaoru Peng;Miaoru Peng;Yuanshu Zhao;Yuanshu Zhao;Shuxian Li;Shuxian Li;Ping Li;Ping Li;Fan Zhang;Fan Zhang;Xiaodong Fu;Xiaodong Fu;Xiaoyang Xu;Xiaoyang Xu et al., Niacin-induced lysosomal free cholesterol efflux via LXRα-mediated signaling pathways in macrophages retards the progression of atherosclerosis, , (2025)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.