

HSF2BP Rabbit Polyclonal Antibody

Catalog No: #53190

Package Size: #53190-1 50ul #53190-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HSF2BP Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of mouse HSF2BP (NP_083178.1).
Conjugates	Unconjugated
Other Names	HSF2BP
Accession No.	Swiss Prot:_Gene ID:74377
Calculated MW	37kDa
SDS-PAGE MW	38kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Western blot analysis of extracts of mouse testis, using HSF2BP Antibody.



Immunofluorescence analysis of L929 cells using HSF2BP Rabbit pAb.

Background

HSF2 binding protein (HSF2BP) associates with HSF2. The interaction occurs between the trimerization domain of HSF2 and the amino terminal hydrophilic region of HSF2BP that comprises two leucine zipper motifs. HSF2BP may therefore be involved in modulating HSF2 activation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.