

Histone H2AX Rabbit Polyclonal Antibody

Catalog No: #53256

Package Size: #53256-1 50ul #53256-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Histone H2AX Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human Histone H2AX.
Conjugates	Unconjugated
Other Names	H2A.X;H2A/X;H2AX;Histone H2AX;H2AFX;histone H2AX;gamma H2A.X;αH2AX;H2AX
Accession No.	Swiss Prot:P16104Gene ID:3014
Calculated MW	15kDa
SDS-PAGE MW	15kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:2000 - 1:5000IHC 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using Histone H2AX antibody.



Immunohistochemistry of paraffin-embedded mouse brain using Histone H2AX antibody.



Immunohistochemistry of paraffin-embedded rat heart using Histone H2AX antibody.



Immunohistochemistry of paraffin-embedded rat brain using Histone H2AX antibody.

Background

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene encodes a replication-independent histone that is a member of the histone H2A family, and generates two transcripts through the use of the conserved stem-loop termination motif, and the polyA addition motif.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.