

RNMT Rabbit Polyclonal Antibody

Catalog No: #53307

Package Size: #53307-1 50ul #53307-2 100ul


Description

Product Name	RNMT Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of mouse RNMT (NP_080716.1).
Conjugates	Unconjugated
Other Names	RNMT;CMT1;CMT1c;MET;Met;RG7MT1;cm1p;hCMT1;hMet;2610002P10Rik;AI848273;Rg7mt1;mKIAA0398
Accession No.	Swiss Prot:Q9D0L8Gene ID:67897
Calculated MW	41kDa/46kDa/53kDa
SDS-PAGE MW	55kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using RNMT antibody.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.