

USP30 Rabbit Polyclonal Antibody

Catalog No: #53436

Package Size: #53436-1 50ul #53436-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	USP30 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human USP30 (NP_116052.2).
Conjugates	Unconjugated
Other Names	USP30
Accession No.	Swiss Prot:Q70CQ3Gene ID:84749
Calculated MW	58kDa
SDS-PAGE MW	70kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Western blot analysis of extracts of MCF-7 cells, using USP30 antibody.

Background

USP30, a member of the ubiquitin-specific protease family (see USP1, MIM 603478), is a novel mitochondrial deubiquitinating (DUB) enzyme (Nakamura and Hirose, 2008 [PubMed 18287522]).[supplied by OMIM, Dec 2008]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.