

# YTHDF1 Rabbit Polyclonal Antibody

Catalog No: #53535

Package Size: #53535-1 50ul #53535-2 100ul


## Description

|                       |                                                           |
|-----------------------|-----------------------------------------------------------|
| Product Name          | YTHDF1 Rabbit Polyclonal Antibody                         |
| Host Species          | Rabbit                                                    |
| Clonality             | Polyclonal                                                |
| Isotype               | IgG                                                       |
| Purification          | Affinity purification                                     |
| Applications          | WB,IF                                                     |
| Species Reactivity    | Human,Mouse,Rat                                           |
| Immunogen Description | Recombinant fusion protein of human YTHDF1 (NP_060268.2). |
| Conjugates            | Unconjugated                                              |
| Other Names           | YTHDF1;C20orf21                                           |
| Accession No.         | Swiss Prot:Q9BYJ9Gene ID:54915                            |
| Calculated MW         | 20kDa/60kDa                                               |
| SDS-PAGE MW           | 61kDa                                                     |
| Formulation           | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.   |
| Storage               | Store at -20°C. Avoid freeze / thaw cycles.               |


## Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


## Images



Western blot analysis of extracts of various cell lines, using YTHDF1 antibody.



Immunofluorescence analysis of L929 cells using YTHDF1 Rabbit pAb.



Immunofluorescence analysis of U2OS cells using YTHDF1 Rabbit pAb.

## Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.