

COL20A1 Rabbit Polyclonal Antibody

Catalog No: #53722

Package Size: #53722-1 50ul #53722-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	COL20A1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human COL20A1 (NP_065933.2).
Conjugates	Unconjugated
Other Names	COL20A1;bA261N11.4
Accession No.	Swiss Prot:Q9P218Gene ID:57642
Calculated MW	135kDa/137kDa/138kDa
SDS-PAGE MW	Refer to figures
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images



Immunofluorescence analysis of rat kidney using COL20A1 antibody.



Immunofluorescence analysis of human placenta using COL20A1 antibody.



Immunofluorescence analysis of mouse brain using COL20A1 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.