

GSDMC Rabbit Polyclonal Antibody

Catalog No: #53877

Package Size: #53877-1 50ul #53877-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GSDMC Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human GSDMC (NP_113603.1).
Conjugates	Unconjugated
Other Names	GSDMC;MLZE
Accession No.	Swiss Prot:Q9BYG8Gene ID:56169
Calculated MW	57kDa
SDS-PAGE MW	58kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Western blot analysis of extracts of HeLa cells, using GSDMC antibody.



Immunofluorescence analysis of L929 cells using GSDMC antibody.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.