

ZNF624 Rabbit Polyclonal Antibody

Catalog No: #53945

Package Size: #53945-1 50ul #53945-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF624 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZNF624 (NP_065838.2).
Conjugates	Unconjugated
Other Names	ZNF624
Accession No.	Swiss Prot:Q9P2J8Gene ID:57547
Calculated MW	85kDa/99kDa
SDS-PAGE MW	—
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.