

CSK Rabbit Polyclonal Antibody

Catalog No: #54300

Package Size: #54300-1 50ul #54300-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CSK Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CSK (NP_001120662.1).
Conjugates	Unconjugated
Other Names	CSK
Accession No.	Swiss Prot:P41240Gene ID:1445
Calculated MW	50kDa
SDS-PAGE MW	50kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using CSK Rabbit pAb.



Immunofluorescence analysis of L929 cells using CSK Rabbit pAb.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.