

CNNM1 Rabbit Polyclonal Antibody

Catalog No: #54377

Package Size: #54377-1 50ul #54377-2 100ul


Description

Product Name	CNNM1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CNNM1 (NP_065081.2).
Conjugates	Unconjugated
Other Names	CNNM1;ACDP1;CLP-1
Accession No.	Swiss Prot:Q9NRU3Gene ID:26507
Calculated MW	—
SDS-PAGE MW	104kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using CNNM1 Rabbit pAb.



Immunofluorescence analysis of C6 cells using CNNM1 antibody.



Immunofluorescence analysis of L929 cells using CNNM1 antibody.



Immunofluorescence analysis of U-2 OS cells using CNNM1 antibody.

Background

This gene encodes a member of the ancient conserved domain protein family. The encoded protein may bind copper. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.