

CYB5R1 Rabbit Polyclonal Antibody

Catalog No: #54385

Package Size: #54385-1 50ul #54385-2 100ul


Description

Product Name	CYB5R1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human CYB5R1 (NP_057327.2).
Conjugates	Unconjugated
Other Names	CYB5R1;B5R.1;B5R1;B5R2;NQO3A2;humb5R2
Accession No.	Swiss Prot:Q9UHQ9Gene ID:51706
Calculated MW	—
SDS-PAGE MW	34kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using CYB5R1 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Human esophageal using CYB5R1 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Rat ovary using CYB5R1 Rabbit pAb.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.