

CAMKV Rabbit Polyclonal Antibody

Catalog No: #54395

Package Size: #54395-1 50ul #54395-2 100ul


Description

Product Name	CAMKV Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CAMKV (NP_076951.2).
Conjugates	Unconjugated
Other Names	CAMKV;1G5;VACAMKL
Accession No.	Swiss Prot:Q8NCB2Gene ID:79012
Calculated MW	51kDa/54kDa
SDS-PAGE MW	54kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using CAMKV Rabbit pAb.



Immunofluorescence analysis of C6 cells using CAMKV antibody.



Immunofluorescence analysis of HeLa cells using CAMKV antibody.



Immunofluorescence analysis of L929 cells using CAMKV antibody.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.