

CAMSAP1 Rabbit Polyclonal Antibody

Catalog No: #54467

Package Size: #54467-1 50ul #54467-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CAMSAP1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human CAMSAP1 (NP_056262.3).
Conjugates	Unconjugated
Other Names	CAMSAP1
Accession No.	Swiss Prot:Q5T5Y3Gene ID:157922
Calculated MW	—
SDS-PAGE MW	240kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of HeLa cells, using CAMSAP1 Rabbit pAb.



Immunofluorescence analysis of L929 cells using CAMSAP1 antibody.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.