

ARHGAP30 Rabbit Polyclonal Antibody

Catalog No: #54468

Package Size: #54468-1 50ul #54468-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARHGAP30 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ARHGAP30 (NP_859071.2).
Conjugates	Unconjugated
Other Names	ARHGAP30
Accession No.	Swiss Prot:Q7Z6I6Gene ID:257106
Calculated MW	—
SDS-PAGE MW	119kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images



Western blot analysis of extracts of Mouse spleen, using ARHGAP30 Rabbit pAb.



Immunofluorescence analysis of HeLa cells using ARHGAP30 Rabbit pAb.



Immunofluorescence analysis of NIH3T3 cells using ARHGAP30 Rabbit pAb.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.